

Springfuser Flow Rate Chart

Product Code	Flow Rate	Flow Rate when using 30ml Syringe & Pump
SPRI T 5	10mL/5 mins	30ml/15 mins
SPRI T 15	10mL/15 mins	30ml/45 mins
SPRI T 30	10mL/30 mins	30ml/90 mins
SPRI T 1	10mL/60 mins	30ml/3 hrs
SPRI T 2	10mL/2 hrs	30ml/6 hrs
SPRI T 4	10mL/4 hrs	30ml/12 hrs
SPRI T 5	10mL/5 hrs	30ml/15 hrs
SPRI T 13	10mL/8 hrs	30ml/24 hrs
SPRI T 12	10mL/12 hrs	30ml/36 hrs
SPRI T 24	10mL/24 hrs	30ml/72 hrs
SPRI T 48	10mL/48 hrs	30ml/144 hrs
SPRI T 72	10mL/72 hrs	30ml/216 hrs

Spring Fuser Pump 10ml - SPRI S 30
 Spring Fuser Pump 30ml - SPRI S 301
 Spring Fuser 30ml Syringe - SPRI S 4

Note: All tubing is supplied with a 10ml syringe and labelled for 10ml flowrate. When using a 30ml Springfuser Pump purchase a 30ml Springfuser Syringe (SPRI S 4) and tubing. Therefore if you use 10/5min tubing you deliver 30ml over 15 mins.

