

KRUUSE Rehab Lifting Harness

The KRUUSE Rehab Lifting Harness for front or hind legs, lifts and supports animals with age-related problems, injuries and mobility issues. Ranging from neurological issues, structural deficiencies, general weakness and the need to safely navigate stairs.

Benefits & Features

- Soft and pliable handles for heightened comfort
- Handles allow free passage of straps for improved mobility
- Touch fastener for safe and adjustable fit
- Reflective tape for better visibility
- Adjustable straps to help the pet owner stand up straight while handling the pet patient

KRUUSE Rehab Lifting Harness is an effective walking aid for animals with

- Arthritis
- Hip dysplasia
- Spinal trauma
- General fatigue

The KRUUSE Rehab Lifting Harness helps prevent back injury to pet owners from lifting pets improperly. The harness also supports cats and dogs when nature calls.




Size Guide

KRUUSE Rehab Lifting Harness for Front Legs

- C Measure the circumference of the chest directly behind the front legs.
- D Measure the circumference of the front leg at the point where it meets the body.


Provet Code	Size	Lifting Harness for Front Legs	
		C	D
		Chest	Front leg
HARN K L01	XS	40-46 cm	24-26 cm
HARN K L02	S	44-52 cm	28-32 cm
HARN K L03	M	52-62 cm	34-38 cm
HARN K L04	L	62-74 cm	44-48 cm
HARN K L05	XL	74-86 cm	54-58 cm



KRUUSE Rehab Lifting Harness for Hind Legs

- A Measure the circumference of the abdomen directly in front of the back legs.
- B Measure the circumference of the upper thigh (back leg) at the point where it meets the body.

Provet Code	Size	Lifting Harness for Hind Legs	
		A	B
		Abdomen	Upper thigh
HARN K L06	XS	30-36 cm	28-32 cm
HARN K L07	S	36-43 cm	33-37 cm
HARN K L08	M	43-51 cm	38-42 cm
HARN K L09	L	51-61 cm	48-52 cm
HARN K L10	XL	61-73 cm	54-58 cm



38% Nylon
62% Neoprene



To order, contact your local Customer Service Team, Provet Business Account Manager or order online through Provet Connect connect.provet.com.au